



# CONSUMABLES FOR SPECTROSCOPIC SAMPLE PREPARATION

SPECTROSCOPY

# SPECTROSCOPIC SAMPLE PREPARATION

Sample preparation of metals and materials have become more and more important because of the rapid development and improvement of both software as well as OES and XRF-devices during the past few years that shifts the detection limit for trace analyses

It is crucial to have the sample properly prepared. The sample needs to be both representative, homogeneous and with an even surface in order to eliminate factors that can influence the results.



## SPECTROSCOPY

## DEMPAX

Paper Grinding Disc for SPECTRAL 250 &amp; 350

Order No	Code	Type	Diameter	Grit Size	Quantity /Pack
37-040-060-22	DEMPAX-S	Corundum	250	60	20
37-065-060-40	DEMPAX-S	Corundum	350	60	20
36-040-060-22	DEMPAX-SZ	Zirconium oxide	250	60	20
36-065-060-40	DEMPAX-SZ	Zirconium oxide	350	60	20

## GRINDING STONES

for SPECTRAL PG 52/AG 102

Order No	Code	Description	Grit Size	Quantity / Pack
80-150	GSW 60	AL <sub>2</sub> O <sub>3</sub> Grinding Stone for medium hard steels	60	1
80-151	GSR 60	AL <sub>2</sub> O <sub>3</sub> Grinding Stone for cast iron and hard steels	60	1
80-155	GSW 36	AL <sub>2</sub> O <sub>3</sub> Grinding Stone for medium hard steels	36	1
80-156	GSR 36	AL <sub>2</sub> O <sub>3</sub> Grinding Stone for cast iron and hard steels	36	1

## TIPO

Milling Tips for SPECTRAL MM 102

Order No	Code	Description	Quantity / Pack
80-204	TIPO	Milling tips for non-ferrous materials	[1 set = 10 pcs.]
80-205	TIPO	Milling tips for cast irons	[1 set = 10 pcs.]
80-206	TIPO	Milling tips for medium hard and hard steels materials	[1 set = 10 pcs.]
80-207	TIPO	Ceramic Milling tips for cast iron & steel materials (GR 1838 & GR 1839)	[1 set=10 pcs.]



┌

Distributor:

└

CMS-TB-01